

lumiNIR

Congratulations on your lumiNIR purchase!

Follow these simple steps for maximum results:

- Read this Quick Start Guide
- Take note of the treatment area before you begin to measure results
- Take a “before” picture if applicable

Treatment steps:

- Cleanse and dry skin
- Turn on device and select intensity level (see chart on back)
- Place light directly on clean, dry skin
- Device will beep when treatment time ends
- Move on to your next treatment area

How to charge your lumiNIR:

The lumiNIR wand uses a standard Micro USB charging cable. This can be combined with a wall plug attachment similar to that used by many mobile devices. When charging your wand, the red light will change to green in approximately two hours. **NOTE:** Overcharging the wand (keeping the wand on the charger past the indicator light turning green), will shorten the lifespan of the wand and require a replacement.

Get to know your device:



lumiNIR Intensity Options

The lumiNIR offers four intensity levels. Level 1 is considered one dose.

This is the minimum amount needed for a treatment. You can increase the dosage based on your need.

Intensity Level	Light Indicator <small>(the light on the on device will change color for each level)</small>	Treatment Time
Level 1 (default)	Green light	35 seconds
Level 2	Blue light	1:10 minutes
Level 3	Green/Red light	1:45 minutes
Level 4	Red light	2:20 minutes

Cleaning:

Use a damp, soft, lint-free cloth and water or rubbing alcohol to clean the LED attachments. Makeup wipes are also a good option.

Do not immerse in water at any time. The back of the LED heads must remain dry.

For additional information, visit [sunlighten.com/customer-service/luminir-support](https://www.sunlighten.com/customer-service/luminir-support).

Scan the QR code with your smart phone to access the above URL.

